



# CFPM300 & CFPM301 Imaging Tables

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures.



Integrated Headrest Table

## Motions

- Motion 1 - Height Range 26" - 44"
- Motion 2 - Trendelenburg Tilt +15/-12°
- Motion 3 - Longitudinal Travel 10"  
(CFPM301)
- Motion 3 - Lateral Tilt  $\pm 15^\circ$  (CFPM300)

## Features

- Choice of integrated headrest or rectangular carbon fiber tops with 56" of metal free imaging
- 2" Radiolucent Comfort Foam<sup>®</sup> table top pad
- Warranty is 3 years parts and 2 years labor

# SPECIFICATIONS

- Large Adjustable Height Range 26" - 44"
- 56" Metal Free Imaging Space
- Cost Effective 3 Movement Table
- 500 lb Patient Weight Capacity
- Warranty: 3 years parts, 2 years labor

CFPM300 #81803 Integrated #81810 Rectangular CFPM301 #81804 Integrated #81811 Rectangular

## INCLUDED FEATURES

Carbon Fiber Top Integrated Headrest 22" (56cm) wide x 84" (214cm) long  
56" Metal free imaging space Rectangular 24" (61cm) wide x 80" (203cm) long

## CFPM300 MOVEMENTS

Motion 1 Electric Height Range 26" - 44" (66cm-111cm)  
(4" casters w/o pad)

Motion 2 Electric Trendelenburg tilt +15°/-12°

Motion 3 Electric Lateral tilt ±15°

## CFPM301 MOVEMENTS

Motion 1 Electric Height Range 26" - 44" (66-111cm)  
(4" Casters w/o pad)

Motion 2 Electric Trendelenburg tilt +15°/-11°

Motion 3 Electric Longitudinal travel 10" (25cm)

Controls All functions hand control

Weight Capacity 500 lb (227kg)

Product Weight 475 lb (216kg)

Included Features 2" (5cm) Radiolucent Comfort Foam™ table top pad, 4" (10cm) locking casters, Crescent face rest pad (included with integrated headrest top), patient restraint straps (2)

Warranty 3 years parts, 2 years labor

Safety Listings FDA listed, CE marked, cETLus listed, CB Certificate

## OPTIONS *must be ordered with table.*

5" (13cm) casters Upgrade from the standard 4" (10cm)

Foot Control All functions control

Non-USA electrical systems Please inquire

## ACCESSORIES *order any time.*

	#
Arm Hammock	Use with Integrated Headrest Top 48751
Carbon Fiber Arm Board (6" wide)	180° Articulation 83575
Fluoro Extender	Polycarbonate 33975
Interventional Supine Pillow System	Contoured pad and support wedge pad 81821
Spine Positioning System II	8 piece system 81819
T-rail Clamp	6" (15cm) long rail 63336



Rectangular Top Table

# ACCESSORIES

## SPINE POSITIONING SYSTEM II

Radiolucency & adjustability for ideal cervical, thoracic & lumbar visualization. Minimize unwanted movement & reduce procedure time.



## ARM HAMMOCK

This unique product will comfortably & easily support your patient's arms when prone. Easily attached to the Integrated Head Rest table top.



## SINGLE CARBON FIBER ARM BOARD

Fully radiolucent with a 180° range of motion, the arm board requires no mounting or attachment hardware and simply slides under the table top pad. It can be used on any fluoroscopy table & will support arms and feet in a variety of positioning configurations.



## FLUORO EXTENDER

Made of radiolucent, polycarbonate material, it's organic shape is conducive to a variety of table placement applications. 30" long x 23.5" wide x 1/4" thick



## 6" ADJUSTABLE T-RAIL CLAMP

The adjustable T-rail clamp can be positioned along the grooved table top edge. Clamp is also removable.



## 5" CASTER (OPTIONAL)

Recommended for mobile service providers.



717.759.3146 [www.oakworksmed.com](http://www.oakworksmed.com) email: [medsales@oakworks.com](mailto:medsales@oakworks.com)